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Product Specifications

Measurements	ECG Dynamic Range	-10mV to +10mV
	Heart Rate (stationary and ambulatory)	30 – 200 Beats per Minute ($\leq \pm 5$ or 10% Beats per Minute, whichever is greater)
	Respiration Rate	<ul style="list-style-type: none"> ○ 10-30 Breaths per Minute with a mean absolute error of less than 3 Breaths per Minute, validated by clinical studies ○ 4-42 Breaths per Minute with a mean absolute error of less than 1.5 Breaths per Minute, validated by simulation studies
	Body Temperature	32°C - 42°C ($\leq 1^\circ\text{C}$, mean absolute error)
	Skin Temperature	15°C – 50°C ($\leq 0.3^\circ\text{C}$)
	Fall Detection	Fall or No Fall (> 90% Sensitivity and >98% Specificity)
	Step Count	< 5% Absolute Error Compared to Manual Count Step count is reset to 0 after step count 65535 is reached.
	Posture Detection	Lying down, Upright, Walking, Running, or Leaning (>70% Accuracy Compared to Visual)
Communications	Bluetooth (BT4.1)	Max. 10 Meters (30 Feet Line of Sight)
	Radio Modulation	FSK (Frequency Shift Keying)
	Radio Frequency	2.4 – 2.5GHz
	Transmit power	$\leq 10\text{dbm}$
	Security	AES-CCM 128 Bit Encryption (Advanced Encryption Standard-CCM mode)
Battery	Battery Type	Zinc Air
	Battery Voltage	DC 1.4V
	Battery Life	120 or 168 Hours
Operating Conditions	Ambient Temperature	10 – 40°C
	Humidity	10 – 95% RH (noncondensing)
	Altitude	<3000 m
	Barometric Pressure	70 kPa to 102 kPa

Electromagnetic Emission Declaration

The VitalPatch device is intended for use in the electromagnetic environment specified below. The end user of the device should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment
RF emissions CISPR 11	Group 1	The VitalPatch device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The VitalPatch device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

FCC Compliance (FCC ID:SPO-VCI-VP2)

The VitalPatch device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation (FCC Title 47, Subpart A, Part 15.19(3)).

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment (FCC Title 47, Subpart A, Part 15.21)

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures (FCC Title 47, Subpart B, Part 15.105(b)):

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Guidance and Declaration – Electromagnetic Immunity

(For ME equipment ME system that are not life-supporting)

The VitalPatch device is intended for use in the electromagnetic environment specified below. The end user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment- guidance
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the VitalPatch device than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:</p> $d = 1.17\sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.33\sqrt{P} \text{ 800MHz to 2.7 GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. **Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people. **Note 3:** UT is the a.c. mains voltage prior to application of the test level. **(a)** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. **(b)** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and VitalPatch

(For ME equipment ME system that are not life-supporting)



The VitalPatch device is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The end user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the VitalPatch device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)	
	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.7 GHz $d = 2.33\sqrt{P}$
0.01	0.17	0.23
0.1	0.37	0.74
1	1.17	2.33
10	3.69	7.38
100	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. The VitalPatch device complies with the applicable requirements and relevant provisions of the Radio Equipment Directive 2014/53/EU (RED). **Note 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. **Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

General Symbols

Symbol	Title	Symbol	Title
IP24	Protected against splashing water	IP27	Protected against submerging in water (up to 1 meter for 30 minutes)
	Re-use is not allowed		Read usage instructions
	Properly dispose of EEE (Electrical and Electronic Equipment)		Non-ionizing radiation

Symbol	Title	Symbol	Title
	Defibrillation proof type CF applied part		MR Unsafe
	Prescription only		Manufacturer
	Caution, consult documents		Not to be used in case package is damaged
	Catalogue number		Batch code
	Use by date		Temperature limits (Storage)
	Humidity limits (Storage)		Contents (Numeral represents quantity of units inside)
	Underwriters Laboratories MEDICAL — PATIENT MONITORING EQUIPMENT AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI ES60601-1 (2005), "Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance; ANSI/AAMI/IEC 60601-2-25, "Medical Electrical Equipment - Part 2-25: Particular Requirements for the Basic Safety and Essential Performance of Electrocardiographs" E358758		

Document History

Revision	Description
Rev. M	Initial release of this document.
Rev. N	Added instructions for performing posture calibration. Revised instructions for body temperature calibration to include a recommended schedule for calibration.
Rev P	Add instruction for 120 or 168 hours of battery life. Minor changes to Step 4 Position VitalPatch on the body for clarity. Update graphics.
Rev R	Add support email to contacts page
Rev S	Remove CE and EC rep symbols in general symbols table. Remove EC rep address from page 13 contact info